



# Proven Solutions for DRIVE-THRU TESTING

## Challenge

Optimize COVID-19 Drive-Thru specimen collection with safety speed and accuracy

## Solution

Leverage mobile barcode technology for data collection, positive patient ID, accurate labeling and printing.

### GOOD

## Laptop-Based/Electronic Health Record Configuration

Aprox. 5-6 minutes | No physical contact | Additional Footsteps | Requires Internet Connection



Patients Line Up



ID scanned patient data captured



Healthcare Technician leaves patient to enter additional data into laptop



Labels printed



Healthcare Technician returns to patient



Swab Taken



Patient Exits

### Better

## Turn-Key Mobile Configuration

Under 3 minutes | No Physical Contact | No Manual Entry | Not internet connection



Patients Line Up



ID scanned patient data captured



Labels Printed



Swab Taken



Patient Exits

### PLUS

## Emergency Wristband Kit

For Patients in Critical Need  
Can Be Used with any Configuration



ID scanned patient data captured



Barcode Wristband ID Printed



Patient immediately taken to hospital care